



INTERNATIONAL PRE-LAUNCH NOTIFICATION

Product(s) Name:	<i>Carriere® Motion Clear™</i> Class II Appliance
Brand:	<i>Carriere</i>
Anticipated Launch Date:	Q2 2016

Overview

We are excited to announce we will be launching our *Carriere Motion Clear* Class II appliance to the international market at the end of June 2016!!! The *Motion Clear* Class II appliance is an elegant, aesthetic, and minimally invasive solution for correcting Class II dental malocclusions in patients of all ages. Like its metal counterpart, the new appliance is used at the beginning of treatment, before braces or aligners, and is designed to establish a Class I occlusion in only 3-4 months...faster than any other appliance on the market today*! The *Motion Clear* Class II appliance is ideal for even the most aesthetically demanding Class II patients, including Invisalign and ceramic bracket cases.



The *Motion Clear* Class II appliance bonds directly, via light-cure adhesive, to the maxillary canine and first or second permanent molar. A solid and consistent source of anchorage on the mandible must be selected to avoid protrusion of the lower incisors. Sources of anchorage are based on clinician preference (ie lingual appliance, mini-screws), however an upper Essix® Appliance** with direct bond tubes on lower molars is recommended. Intra-oral elastics activate the maxillary posterior segment and move it into an ideal Class I position. Once the Class I occlusion is achieved, the clinician is able to rapidly finish the case using braces or aligners.

Clinical Applications:

Treating a full range of class II malocclusions, including:

- Division 1 & Division 2 Class II cases (Symmetrical & Asymmetrical)
- Phase I of mixed Class II patients with erupted molars
- Class I or II extraction cases
- Unilateral Class II cases
- Class I cases with mesially positioned maxillary molars



Patented Pad Design

Key Product Features:

- The aesthetic clear arm is comprised of an advanced medical grade polymer for excellent strength and durability. The molar pad is made of stainless steel.
- Resistant to staining, including coffee, red wine, and soft drinks
- Patented anterior pad with special dovetailed grooves provides optimal retention and bond strength
- The ball and socket joint provides maximum freedom of movement with built-in stops to prevent unwanted over-rotation and/or tipping
- Available in six (6) sizes (16 mm, 18 mm, 20 mm, 23 mm, 25 mm, 27 mm) and color-coded for easy identification



Key Product Benefits

The *Motion Clear* Class II appliance provides clinicians and patients access to a truly complete aesthetic solution for correcting Class II malocclusions, resulting in:

- Shorter treatment times
- Happier patients
- Greater referrals
- Great revenue and higher profit margins

Part Numbers:

The *Motion Clear* Class II appliance is available in Left and Right configurations, and offered in singles, sets of two, a Trial kit (3 sets), an Intro kit (6 sets), and a Standard kit (20 sets). Pricing is TBD and will be communicated prior to launch. Corresponding part numbers are listed below.

DESCRIPTION/SIZES	COLOR	LEFT	RIGHT	SET
Carriere Motion Clear Class II 16 mm	○	424-816LC	424-816RC	424-816C
Carriere Motion Clear Class II 18 mm	●	424-818LC	424-818RC	424-818C
Carriere Motion Clear Class II 20 mm	●	424-820LC	424-820RC	424-820C
Carriere Motion Clear Class II 23 mm	●	424-823LC	424-823RC	424-823C
Carriere Motion Clear Class II 25 mm	●	424-825LC	424-825RC	424-825C
Carriere Motion Clear Class II 27 mm	●	424-827LC	424-827RC	424-827C
Carriere Motion Clear Class II Trial Kit (22, 25, 27mm) 1-set each				424-800C
Carriere Motion Clear Class II Intro Kit (6 sets)				424-801C
Carriere Motion Clear Class II Standard Kit (20 sets)				424-802C

Product Registration

All non-European & non-Canadian distributors intending to carry the *Motion Clear* Class II appliance each must register the products for sale according to their local rules and regulations prior to receipt of product from *Ortho Organizers*®. Each distributor must also submit a copy of their local registration certificate to *Ortho Organizers* before their initial order can be released and shipped to ensure expeditious delivery.

Starting the registration process in advance will greatly improve the efficiency and timeliness of product availability in each international location. Please contact Colleen Boswell- Director of Regulatory & Quality Affairs (Colleen.Boswell@henryschein.com) to obtain the necessary information to register the products listed in this notification.

* "Treatment effects of the *Carriere* Distalizer (*Motion*) using lingual arch and full fixed appliances", *Journal of the World Federation of Orthodontists*- May 2014

**Dentsply Raintree Essix, Sarasota, FL, USA